## DOASENSE becoming a standard of care in Italy and Estonia

### DOASENSE is the recommended Point Of Care Test for the prescription of Andexanet alpha (ONDEXXYA) in Tuscany (Italy)

In August 2023 the highest regional health authorities in Tuscany (Italy) enacted a new procedure for the appropriate prescription of ONDEXXYA that is recommended to all professionals involved. It clearly stipulates DOASENSE and provides detailed guidance.

Document 1: The letter sent to all Emergency Departments stipulates that: “The procedure provides that, in the decision-making process for the prescription of the ONDEXXYA drug, the use of the following test is indicated: test strip for DOACs – DOASENSE Dipstick Ref 001

Further, the letter provides this recommendation: “We hereby recommend that the companies and organisations in question make available the diagnostic indicated above in all the Emergency Urgent Structures of the Region, together with the specific reader.”

Document 2: This letter was sent together with an order form that clinicians must use to prescribe ONDEXXYA. The section “DOSAGE NAO” (Dose if DOAC) lists the 2 options for testing, i.e. plasma level >30ng/mL or Urinary (POCT) positive.

Document 3 was also attached to the communication on indications for use of Andexanet Alfa. Page 8 provides an algorithm to support the prescription of the drug showing the same 2 alternatives for the testing: plasma level >30ng/mL or Urinary (POCT) positive.

While Andexanet Alpha is not yet reimbursed in Italy, its prescription is now conditional to the adherence to this protocol in Tuscany. This is an important milestone towards the adoption of DOASENSE as a standard of care to determine the presence of clinically relevant level of DOACs in order to support clinicians with additional information otherwise not available in emergency settings.

Attached are the original documents in Italian and their non-certified translations, as obtained through a smart translator.

### DOASENSE is also now imbedded in anaesthesia guidelines in Estonia

In the guidelines for "Anticoagulation drugs management in elective and emergency surgery" from the Estonian Society of Anesthesiology. (<https://anest.ee/ravijuhised/#> then click on : "Uus! Antitrombootilised ravimid plaanilises ja erakorralises kirurgias")

DOASENSE is listed on Table 10 (page 14): as a clear option to exclude the presence of Dabigatran, Rivaroxaban, Apixaban and Edoxaban. The original document and its non-certified translation are attached as Document 4.

Those are new examples of the momentum gained establishing DOASENSE as a standard of care.

**A.I.**

Document set from the regional health authorities in Tuscany, Italy

*ENGLISH Convenience Translation*

**TUSCANY REGION**  
**Regional Council**

**DIRECTION**

**HEALTH, WELFARE AND SOCIAL COHESION**

**PHARMACEUTICAL ASSISTANCE AND DEVICES SECTOR**

*Responsible manager:* **Claudio Marinai**  
claudio.marinai@regione.toscana.it

---

**Subject:** Interruption of availability of the medicinal specialty ONDEXXYA.

General Directorates of the  
Local Health Authorities and AOU of Tuscany  
**ESTAR General Director**  
Director of the Monasterio Foundation

and PC      Director of the Health, Welfare and Social  
Cohesion  
Directorate Dr. Federico Gelli

Head of the Hospital Care, Quality and Clinical  
Networks Sector Dr. \_\_\_\_\_  
Michela Maielli

PEC

Following the previous Circular of this Sector (AOOGRT/PD prot. n. 0250545 of 05/31/2023) with which the availability of the drug ONDEXXYA was ceased, in agreement with the OTGC an uncovered clinical need was identified adequately by the therapies proposed by Aifa with Circular 20012/2023 and already available (PCC 4 palm factors).

To guarantee the citizens of Tuscany all necessary therapies, even if these include the use of a drug, as mentioned in the aforementioned Circular, which is not currently paid for by the NHS, the undersigned Sector and the OTGC have developed the attached procedure. This procedure limits the use of the ONDEXXYA drug only to cases not covered by already available therapies.

Compliance with the indications of the attached procedure is of fundamental importance for the appropriate prescription of the ONDEXXYA drug and therefore compliance is recommended to all professionals involved.

The procedure provides that in the decision-making process for the use of the ONDEXXYA drug the use of the following test is indicated: Estar code:

50395392 Description: test strip for doacs – doasense dipstick REF 0001

We hereby recommend that the companies and organizations in question make available the diagnostic indicated above in all the Emergency Urgent Structures of the Region, together with the specific reader.

Best regards.

The OTGC Director  
Dr. Stefano Grifoni

The Responsible Manager  
Dr. Claudio Marinai

Regione Toscana

## Request for andexanet alfa (Ondexxya®)

All

## for reversal of anticoagulant therapy in case of potentially fatal or uncontrolled bleeding

The undersigned Dr. \_\_\_\_\_ telephone number \_\_\_\_\_

SOD \_\_\_\_\_ Cost Center \_\_\_\_\_

## requests the supply of the drug andexanet alfa (Ondexxya®) 200 mg

indicated for adult patients treated with apixaban or rivaroxaban, when reversal of anticoagulant therapy is required due to life-threatening or uncontrolled bleeding

for the patient (name and surname) \_\_\_\_\_ date of birth \_\_\_\_\_

| Inclusion and exclusion criteria for treatment                                                                                       |                                                                                                                                                                                                                                                                                                                                                            |                                               |                                   |
|--------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|-----------------------------------|
| o MANIFEST ACTIVE bleeding in patient receiving therapy with APIXABAN or RIVAROXABAN<br>o NOT for pre-treatment in emergency surgery |                                                                                                                                                                                                                                                                                                                                                            |                                               |                                   |
| Select at least one condition or                                                                                                     |                                                                                                                                                                                                                                                                                                                                                            |                                               |                                   |
| A ÿ Critical site bleeding (select option)                                                                                           | CEREBRAL: intraparenchymal; subdural; epidural; subarachnoid*<br>o Others affecting the CNS: intraocular; intra or extraaxial spinal<br>o HEMOPERICARDIUM with signs of cardiac tamponade<br>o AIRWAY (including posterior epistaxis), HEMOTHORAX^<br>or INTRAABDOMINAL, RETROPERITONEAL^<br>*GCS>7 and/or blood volume >60 ml<br>^only IF WITH B and/or C |                                               |                                   |
|                                                                                                                                      |                                                                                                                                                                                                                                                                                                                                                            |                                               |                                   |
| B ÿ Hemodynamic instability (select option)                                                                                          | o Systolic BP < 90 mmHg or reduction in Systolic BP values > 40 mmHg<br>o Shock index (HR/BP) > 0.9<br>o Lactates > 4 mmol/L<br>o Urine < 0.5 ml/Kg/min                                                                                                                                                                                                    |                                               |                                   |
| C ÿ Hb reduction ÿ2 g/dl or transfusion of ÿ2 units                                                                                  | If known, specify bleeding site:                                                                                                                                                                                                                                                                                                                           |                                               |                                   |
|                                                                                                                                      |                                                                                                                                                                                                                                                                                                                                                            |                                               |                                   |
| DOSAGE NAO<br>(select the option and attach reference to this document)                                                              |                                                                                                                                                                                                                                                                                                                                                            | o Plasma: >30 ng/ml                           |                                   |
|                                                                                                                                      |                                                                                                                                                                                                                                                                                                                                                            | o Urinary (POCT): POSITIVE                    |                                   |
| Doses for apixaban/rivaroxaban reversal                                                                                              |                                                                                                                                                                                                                                                                                                                                                            |                                               |                                   |
| Inhibitor of FXa                                                                                                                     | Last dose                                                                                                                                                                                                                                                                                                                                                  | Time since last dose before starting Ondexxya |                                   |
|                                                                                                                                      |                                                                                                                                                                                                                                                                                                                                                            | < 8 hours                                     | ÿ 8 hours                         |
| apixaban/<br>rivaroxaban                                                                                                             | ÿ 5 mg apixaban/10 mg rivaroxaban Low dose (5 vials)                                                                                                                                                                                                                                                                                                       |                                               | Low dose (5 vials)                |
|                                                                                                                                      | > 5 mg apixaban/10 mg rivaroxaban High dose (9 vials)                                                                                                                                                                                                                                                                                                      |                                               |                                   |
| Prescribed dosage                                                                                                                    |                                                                                                                                                                                                                                                                                                                                                            |                                               |                                   |
| No. of 200 mg vials                                                                                                                  | Initial intravenous bolus                                                                                                                                                                                                                                                                                                                                  |                                               | IV infusion                       |
| or 5                                                                                                                                 | 400 mg at a target rate of 30 mg/min                                                                                                                                                                                                                                                                                                                       |                                               | 4 mg/min for 120 minutes (480 mg) |
| or 9                                                                                                                                 | 800 mg at a target rate of 30 mg/min                                                                                                                                                                                                                                                                                                                       |                                               | 8 mg/min for 120 minutes (960 mg) |

Regione Toscana

**Request for andexanet alfa (Ondexxya®)**

**for the reversal of anticoagulant therapy in case of  
potentially fatal or uncontrolled bleeding**

All \_\_\_\_\_

Date \_\_\_\_\_

STAMP AND SIGNATURE OF THE REQUESTING DOCTOR

## Indications for the use of Andexanet alfa

### SUMMARY

|                                            |    |
|--------------------------------------------|----|
| 1 PURPOSE.....                             | 2  |
| 2 FIELD OF APPLICATION.....                | 2  |
| 3 REFERENCES.....                          | 2  |
| 4 MECHANISM OF ACTION.....                 | 2  |
| 5 THERAPEUTIC INDICATIONS.....             | 3  |
| 6 DOSAGE and METHOD of ADMINISTRATION..... | 8  |
| 7 HANDLING INSTRUCTIONS.....               | 9  |
| 8 SIDE EFFECTS.....                        | 10 |

## 1 PURPOSE

The purpose of the following work is to outline a guide to using the drug andexanet alfa (Ondexxya®) in the Tuscany Region. Andexanet alfa (Ondexxya®), powder for solution for infusion, is a modified recombinant derivative of factor recombinant DNA technology in Chinese Hamster Ovary cells (Chinese Hamster Ovary, CHO).

## 2 FIELD OF APPLICATION

The procedure is applied within the Emergency Department to achieve an effect *reverse* in patients treated with apixaban and rivaroxaban who experience bleeding potentially fatal or critically ill (see figure 1).

## 3 REFERENCES

- *CPR*
- *AIFA determination of classification of the drug for reimbursement purposes. AIFA determination n. 20012/2023 inclusion of PROPLEX in the list of medicines payable entirely by the NHS*
- *Important information note agreed with the European regulatory authorities and the Italian Medicines Agency (AIFA)*
- *Steffel J, Collins R, Antz M, et al. 2021 European Heart Rhythm Association Practical Guide on the Use of Non-Vitamin K Antagonist Oral Anticoagulants in Patients with Atrial Fibrillation. Europace. 2021 Oct 9;23(10):1612-1676. 10.1093/europace/doi: euab065. Erratum in: Europace. 2021 Jun 28;: PMID: 33895845.*
- *Tomaselli G, Mahaffey K, Cuker A, et al. 2020 ACC Expert Consensus Decision Pathway on Management of Bleeding in Patients on Oral Anticoagulants. J Am Coll Cardiol. 2020 Aug, 76 (5) 594–622.*
- *Rossaint et al. Critical Care (2023) 27:80. The European guideline on management of major bleeding and coagulopathy following trauma: sixth edition*
- *Cao D. et al. ESC. Reversal and removal of oral antithrombotic drugs in patients with active or perceived imminent bleeding. European Heart Journal (2023) 00, 1–15*

## 4 MECHANISM OF ACTION

Andexanet alfa is a recombinant form of the protein called human FXa, which was modified and deprived of the enzymatic activity of FXa. The active site serine has been replaced with alanine, to make the molecule incapable of cleaving and activating prothrombin, and the domain of gamma-carboxyglutamic acid (Gla) has been removed, to deprive the protein of its capacity of assembly into the prothrombinase complex, thus eliminating any effect anticoagulant.

Andexanet alfa is a reversal agent specific for FXa inhibitors. The main mechanism of action is the binding and sequestration of the FXa inhibitor. Furthermore, it was observed that andexanet

alpha binds to TFPI and inhibits it. Inhibition of TFPI activity may increase tissue factor-initiated thrombin generation, inducing a procoagulant effect.

The effects of andexanet alfa can be measured by pharmacodynamic markers, including the free fraction of available FXa inhibitor, and by restoration of thrombin generation.

Commercially available anti-FXa assays are not suitable for measuring anti-FXa activity after administration of andexanet alfa. Due to the reversible binding of andexanet alfa to the FXa inhibitor, the high sample dilution currently used in these assays results in dissociation of the inhibitor from andexanet alfa, resulting in the detection of falsely high anti-FXa activity levels and resulting in therefore a substantial underestimate of the reversal activity of andexanet alfa. In studies determining the

Prospective, randomized, placebo-controlled trials conducted in healthy subjects determined the dose and dosing regimen of andexanet alfa required to reverse anti-FXa activity and restore thrombin generation for FXa inhibitors (apixaban or rivaroxaban ) using modified assays not available on the market.

Maximum reversal of anti-FXa activity was achieved within two minutes of completion of bolus administration. Administration of andexanet alfa as a bolus followed by continuous infusion resulted in a sustained reduction in anti-FXa activity. Anti-FXa activity returned to placebo levels and above approximately two hours after the end of the bolus or infusion, depending on the dosage. When andexanet alfa was administered as a bolus followed by a continuous infusion, the maximum reduction in unbound FXa inhibitors was rapid (within two minutes of the end of the bolus), was maintained throughout the infusion, and then gradually increased with time, reaching a maximum approximately two hours after the end of the infusion. The restoration of thrombin generation after administration was dependent on the dose and dosing regimen and not related to anti-FXa activity beyond approximately four hours (see section "restoration of thrombin generation"). Plasma TFPI activity was shown to be completely inhibited from 2 minutes to 14.5 hours after bolus administration of andexanet alfa and returned to baseline within 3 days. Tissue factor (TF)-initiated thrombin generation rose immediately above baseline (before anticoagulation) and remained elevated for >20 hours, compared to placebo. The plausibility of a pro-coagulant effect of TFPI inhibition is supported by consecutive and sustained increases in D-dimers, TAT and F1+2.

Studies of andexanet alfa in the presence of direct FXa inhibitors in healthy subjects demonstrated dose-proportional pharmacokinetics in the expected therapeutic dose range for both C<sub>max</sub> and area under the curve (AUC). FXa inhibitors did not alter the pharmacokinetics of andexanet alfa at therapeutic levels.

## **5 THERAPEUTIC INDICATIONS**

For adult patients treated with a direct factor Xa (FXa) inhibitor ( apixaban or **rivaroxaban**), when reversal of anticoagulant therapy is required due to bleeding potentially fatal or uncontrolled.

**However, the non-reimbursability of the drug by the NHS and the literature data a provision at the time of drafting this protocol require its use in one extremely selected case history.**

The possibility of using andexanet alfa must be evaluated following an algorithm that provides three steps (see figure 1) in sequence:

### **1. Inclusion and exclusion criteria**

The administration of Andexanet-alfa is **ONLY** possible for patients undergoing therapy oral anticoagulant with APIXABAN or RIVAROXABAN. It is not prescribable for patients in Edoxaban therapy. Treatment of bleeding emergencies in patients on edoxaban provides, as the only therapeutic option, the administration of PCCs (concentrates of the complex prothrombin) at a dosage of 50 U/kg. Furthermore, the use of andexanet -alfa is NOT possible (e.g current lack of data in this clinical setting) for patients treated with factor inhibitors Xa who need to undergo emergency surgery.

### **2. Type of Bleeding**

Patients with obvious active bleeding and at least one of the three are considered eligible clinical conditions listed below:

A) patients with ~~bleeding in critical locations~~ CEREBRAL (intraparenchymal, subdural; epidural; subarachnoid; intraocular; spinal intra or extraaxial; HEMOPERICARDIUM with signs of cardiac tamponade; ways AIR AND HEMOTHORAX; INTRAABDOMINAL, RETROPERITONEAL.

- at the CEREBRAL level if GCS>7 and/or blood volume<60 ml

- at the AIRWAY and HEMOTORAX level: ONLY if point B or C is also present

Intraluminal gastrointestinal bleeding is NOT considered critical or at risk life, unless there are signs of haemodynamic instability.

A) Patients with **hemodynamic instability**. An increased heart rate may be the first sign of hemodynamic instability due to blood loss. Additionally, a systolic blood pressure (BP) <90 mmHg, a decrease in systolic blood pressure >40 mmHg, or changes in orthostatic blood pressure (systolic blood pressure drop ≥20 mmHg or diastolic blood pressure drop ≥10 mmHg when standing) may indicate hemodynamic instability. The shock index (HR/systolic BP) >0.9 and the pressure

arterial differential  $<40/30$  mmHg are reliable indices of the degree of hemodynamic instability and the need for blood transfusions, as well as surrogate markers for organ perfusion, such as the level of blood lactate, the ultrasound status of filling of the vena cava, the urine output  $<0.5$  ml/kg/h). When practical, continuous invasive measurement of mean arterial pressure is considered superior to peripheral measurement, and a value  $<65$  mmHg serves as a cut-off for hemodynamic instability.

B) **Overt and active bleeding with drop in hemoglobin  $\geq 2$  g/dl or administration of  $\geq 2$  units of packed red blood cells.** Bleeding events causing a drop in hemoglobin  $\geq 2$  g/dL or requiring transfusion of  $\geq 2$  units of red blood cells were associated with a significantly increased risk of mortality. Hematuria and posterior epistaxis associated with anemia requiring blood transfusion may require an initial treatment of mechanical and/or endoscopic hemostasis, if there is no hemodynamic instability.

### 3. Apixaban/Rivaroxaban dosage

As reported in the paragraph "Mechanism of action", andexanet-alfa exerts an action of antidote since it binds to Fxa to replace the circulating anticoagulant drug. Because of this, the possible absence of anticoagulant activity of the drug (defined as antiXa activity under the threshold of 30 ng/ml) does not make the action of the antidote reasonable (at least in a context of necessary selection of case studies, and in the absence of current evidence to the contrary). Therefore the verification of this activity is mandatory for the administration of the drug and in every case must accompany the prescription itself.

The reference method for determining the antiXa activity of the anticoagulant drug is represented by the plasma dosage.

However, the plasma dosage:

- Requires execution times that are often not compatible with the emergency/urgency regime
- It is not available in the laboratories of all hospital facilities
- Even when available, not all laboratories offer 24-hour testing

For this reason, it is useful to consider – in emergency/urgent contexts where a response is necessary as timely as possible – the urinary POCT dosage which offers the possibility of knowing in a few minutes: - type of DOAC taken; - a semi-quantitative response with POSITIVE value or NEGATIVE equivalent to a concentration  $\geq$  or  $< 30$  ng/ml.

The doctor treating the patient with a bleeding event during DOAC must evaluate:

- the anticoagulant drug taken with the usual prescribed dosage
- the time of the last known intake

- or the inability to obtain information on the dosage and timing of intake
- any concomitant therapies with an impact on the risk of bleeding: NSAIDs; antiaggregants; 'over-the-counter' drugs that do not require a prescription
- blood sampling for: complete blood count; renal and liver function; PT, aPTT and fibrinogen (hemorrhage profile prepared in the FA web application).

If there is any doubt about not taking it, or there are uncertainties about the type of drug taken, it is the use of DOAC plasma dosage is indicated, if this is available. Alternatively, it is possible to use point-of-care semi-quantitative measurement devices on urine. If the intake of the DOAC dates back to less than 2-4 hours, the administration of 50 g of activated charcoal orally (unless there are contraindications and if the patient is clinically stable).

In all cases where the algorithm is interrupted, the administration of PCCs is indicated (prothrombin complex concentrates) at 4 factors at a dosage of 50 U/kg.

**In fact, it should be remembered that AIFA has authorized the reimbursement by the NHS of the concentrates of the 4-factor prothrombin complex (Proplex) with the same indications as Andexanet-alfa.**

**PCC4 THEREFORE REPRESENT AN ALTERNATIVE TO ANDEXANET ALFA FOR THE FOLLOWING INDICATIONS:**


- **treatment of adult** patients treated with oral anticoagulants direct inhibitors of factor
- **surgical interventions or invasive maneuvers with a high risk of bleeding to be performed urgently, with timing not compatible with the sole suspension of the anticoagulant; life-threatening or uncontrolled bleeding.**

**USE OF ANDEXANET ALFA IN ASSOCIATION WITH OTHER SUPPORT MEASURES**

The safety of andexanet alfa has not been evaluated in patients receiving prothrombin complex concentrates, recombinant factor VIIa, or whole blood within seven days prior to the bleeding event, because they were excluded from the clinical trials. Treatment with procoagulant factors (e.g., 3- or 4-factor prothrombin complex concentrate (PCC)/activated PCC, recombinant factor VIIa, fresh frozen plasma) and whole blood should be avoided unless absolutely necessary, due to lack of data in association with these treatments.

**After each use, the Pharmacy will request the patient's documentation for the clinical evaluation of the results of the therapy.**

**FIGURE 1**

**6 DOSAGE AND METHOD OF ADMINISTRATION** Each vial

contains 200mg of andexanet alfa; after reconstitution, each mL of solution contains 10 mg of andexanet alfa.

Andexanet alfa is administered as an intravenous bolus at a target rate of approximately 30 mg/min over 15 minutes (low dose) or 30 minutes (high dose), followed by a continuous infusion of 4 mg/min (low dose) or 8 mg/min (high dose) for 120 minutes.

|                   | Initial intravenous bolus            | Continuous intravenous infusion   | Total number of vials from 200 mg necessary |
|-------------------|--------------------------------------|-----------------------------------|---------------------------------------------|
| Dose low          | 400 mg at a target rate of 30 mg/min | 4 mg/min for 120 minutes (480 mg) | 5                                           |
| Dose elevation to | 800 mg at a target rate of 30 mg/min | 8 mg/min for 120 minutes (960 mg) | 9                                           |

**APIXABAN REVERSE**

The recommended dosing regimen of andexanet alfa is based on the dose of apixaban administered to the patient at the time of reversal of anticoagulant therapy and the time since the last dose of apixaban.

| Inhibitor of the FXa | Last dose        | Time since last dose before starting Ondexxya |           |
|----------------------|------------------|-----------------------------------------------|-----------|
|                      |                  | < 8 hours or unknown                          | ≥ 8 hours |
| Apixaban             | ≤ 5 mg           | Low dose                                      | Low dose  |
|                      | > 5 mg/not known | High dose                                     |           |

**RIVAROXABAN REVERSE**

The recommended dosing regimen of andexanet alfa is based on the dose of rivaroxaban administered to the patient at the time of reversal of anticoagulant therapy and the time since the last dose of rivaroxaban.

| Inhibitor of the FXa | Last dose | Time since last dose before starting Ondexxya |           |
|----------------------|-----------|-----------------------------------------------|-----------|
|                      |           | < 8 hours or unknown                          | ≥ 8 hours |

|             |                  |           |          |
|-------------|------------------|-----------|----------|
| Rivaroxaban | ≤ 10 mg          | Low dose  | Low dose |
|             | >10 mg/not known | High dose |          |

## INTERACTION WITH HEPARIN

The use of andexanet alfa before heparinization, e.g. during surgery, you must be avoided, as andexanet induces heparin refractoriness. The use of andexanet as an antidote for heparin or low molecular weight heparin has not been evaluated and is not recommended.

## 7 HANDLING INSTRUCTIONS

Each vial must be reconstituted according to the following instructions: 1.

- Remove the snap cap from each vial.
2. Clean the rubber stopper of each vial with an alcohol swab.
3. Withdraw 20 mL of water for injections using a 20 mL syringe (or more) and a 20 gauge (or larger) needle.
4. Insert the syringe needle into the center of the rubber stopper.
5. Slowly depress the plunger and inject the 20 mL of water for injections into the vial, directing the jet towards the internal wall of the vial to minimize the foam formation.
6. Gently swirl each vial until the powder is completely dissolved. NOT SHAKE the vials to avoid foaming. The dissolution time is approximately three to five minutes per vial.
7. Before administration, visually inspect the reconstituted solution, e.g. exclude the presence of corpuscular material and/or color alterations. Do not use in presence of opaque particles or color changes.
8. For effective reconstitution of the necessary dose and to minimize errors, inject 20 mL of water for injections into each vial before proceeding to the next phase.
9. If stored at room temperature, use within eight hours of reconstitution.

### Administration via timed syringe\_\_\_\_\_

1. After reconstituting all the necessary vials, withdraw the solution from each one using the large syringe (50 mL or larger) equipped with a 20 gauge (or larger) needle.
2. Prepare the bolus and infusion in separate large syringes.
3. Due to the larger volume, bolus and high-dose infusion should be divided into additional syringes (two syringes each for bolus and infusion).

4. To prevent accidental transfer of air, keep the syringe needle pointed up and do not put the syringe down between one sample and the next.
5. Before administration, connect accessories (e.g. extension tubing, in-line filter 0.2 or 0.22 micron polyethersulfone (PES) or equivalent low-bond filter protein, timed syringe).
6. Administer the reconstituted solution at the appropriate rate.
7. Dispose of all used syringes, needles and vials as well as any unused reconstituted solution.

#### Administration via intravenous bags

1. After reconstituting all necessary vials, withdraw the reconstituted solution from each using the large syringe (50 mL or more) equipped with a 20 gauge needle (or bigger).
2. Transfer the reconstituted solution from the syringe into an appropriate IV bag
3. Repeat steps 1 and 2 to transfer the entire bolus and infusion volume into bags iv in PO or PVC. 4. It is recommended to divide the bolus and infusion into two bags separated to ensure correct administration speed. You can also use a single PO or PVC IV bag for the bolus and infusion, but speed must be guaranteed correct infusion rate when switching from bolus to infusion.
4. Before administration, connect accessories (e.g. extension tubing, in-line filter 0.2 or 0.22 micron polyethersulfone (PES) or equivalent low-bond filter protein, IV pump).
5. Administer the reconstituted solution at the appropriate rate. Disposal All syringes, used needles and vials as well as unused residues of the reconstituted solution must be disposed of in accordance with local regulations in force.

Unreconstituted medicine should be stored in a refrigerator (2°C – 8°C).

Andexanet alfa does not necessarily need to be brought to room temperature before reconstitution or before administration to the patient. Adopt aseptic techniques during processing reconstitution procedure.

Chemical and physical stability of the reconstituted medicinal product has been demonstrated for 16 hours at 2°C – 8°C in the immediate packaging vial. If necessary, the reconstituted solution transferred into the IV bag can be stored for another eight hours at room temperature. From the point of view microbiological, the reconstituted product should be used immediately. If it is not used immediately, storage times during use and storage conditions before of use are the responsibility of the user.

## 8 SIDE EFFECTS

Safety was evaluated in clinical studies including 247 healthy subjects who received it administered an FXa inhibitor and 352 patients included in a phase IIIb/IV study (ANNEXA-4) suffering from major acute bleeding and being treated with an FXa inhibitor (mainly apixaban and rivaroxaban). In clinical studies with healthy subjects who were administered a FXa inhibitor and subsequently andexanet alfa, no serious adverse reactions were reported or severe.

In pivotal studies 14-503 and 14-504, treatment with andexanet alfa resulted in statistically significant increase in thrombin generation in healthy subjects subjected to anticoagulant therapy with apixaban or rivaroxaban vs. placebo. In study 14-505 (ANNEXA-4), single arm conducted on patients suffering from major acute bleeding and being treated with a FXa inhibitor (apixaban and rivaroxaban), 11% of patients treated with andexanet alfa had experienced one or more thromboembolic events. The median time to first thromboembolic event is ten day state; 38% of patients with thromboembolic events experienced the mole event first three days.

In fact, it is known that patients treated with an FXa inhibitor are affected by underlying pathologies that predispose them to thrombotic events and reversal of therapy with an FXa inhibitor exposes patients to the thrombotic risk linked to the underlying disease.

Furthermore, an independent pro-coagulant effect of andexanet alfa was demonstrated in studies, mediated by inhibition of tissue factor pathway inhibitor (TFPI), which may constitute a risk for the development of thrombosis. The duration of this effect moles haemorrhagic patients is unknown.

**B.I.**

Document set from the regional health authorities in Tuscany, Italy

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