



DOAC Dipstick



DIAGNOSTIC TEST FOR DOACs IN URINE



EN

Intended Use:

The diagnostic test strip DOAC Dipstick is intended for qualitative detection of the absence or presence of direct oral anticoagulants (DOACs: Dabigatran, Apixaban, Edoxaban and Rivaroxaban) in human urine by visual identification of colours. The DOAC Dipstick is an in vitro diagnostic test intended for professional use only.

Summary and Explanation

The clinical importance of positive results of a DOAC in human urine relates to the presence of DOACs in blood. DOACs are excreted rapidly into urine starting 1 to 2 hours after intake of the medication. A specific and rapid detection indication by a point-of-care test may support diagnosis of anticoagulant therapy with DOACs especially in emergency medicine. Medical decision making may be accelerated. Typical indications are patients with ischemic or haemorrhagic stroke with indication for fibrinolytic therapy or administration of a specific antidote, major traumas, emergency procedures, spontaneous thrombotic and bleeding events during oral anticoagulant therapy, and situations without available medication history. The kidney excretes creatinine and DOACs into urine. Their excretion decreases with impairment of renal function. Consequently, DOACs accumulate in blood with decrease of renal function. Therefore, the parameter of creatinine was included on the DOAC Dipstick. If creatinine is low in urine, false negative results of DOACs may be detected in urine samples using DOAC Dipsticks.

Material Provided With the Test

Test strips: 12
Test strip tube container with printed colour scale and cap: 1
Instructions for use: 1


Materials Required but not Provided

Clean container made of polypropylene for urine sample collection.
Timer.

Principle of the DOAC Test Strip

The test consists of a change of colour upon reaction of Factor Xa or Thrombin on the release of the chromophore bound to a substrate. The release of the chromophore is negatively related to the amount of DOACs in urine. Colours for oral direct Factor Xa and Thrombin inhibitors are different. The colour of the pads on the test strips changes within 10 min and can be identified by naked eye. The colours allow the detection of DOACs in a urine sample with interpretation as "negative" in the absence of a DOAC and as "positive" in the presence of a DOAC. Respective colours for comparison are printed on the test tube containing the test strips.

The DOAC Dipstick has four different test pads used for analysis as follows:



Schematic of the DOAC Dipstick

Principle of Tests

Thrombin Inhibitor – Thrombin releases a chromophore from a Thrombin-specific peptide that is inhibited in the presence of a direct oral Thrombin inhibitor.

Factor Xa Inhibitor – Factor Xa releases a chromophore from a Factor Xa-specific peptide that is inhibited in the presence of a direct oral Factor Xa inhibitor.

Urine Colour – This pad does not contain any reagents and is used for assessing the impact of the colour of the patient urine that was used.

Creatinine – The test is based on reaction of creatinine with 3,5-dinitrobenzoic acid in alkaline medium (Benedict-Behre reaction).

Used test strips must be treated as hazardous waste according to national biohazard and safety guidelines or regulations. All reagents of this kit have been found to be uninfected. However, materials contaminated with human urine can be infectious. For this reason, used test strips should be treated as potential biohazards in use and for disposal. If contamination of clothing occurs: Rinse skin with water or shower. Avoid contact with skin and eye. If skin irritation occurs: Consult a physician in all serious cases of health damage. In case of an accidental ingestion wash out the mouth and drink approximately 0.5 l of water. In case of eye contact, rinse the eye quickly and thoroughly using a stream of clean water.

Procedure Notes Before Performing the Test

Do not use test strips with an expiry date that has already passed. Carefully read the instructions for use before starting the test. The instructions must be followed exactly to obtain accurate results. This test is for professional in vitro diagnostic use only.

Do not touch test pads of the strip. Handle the test strips only at the handling area opposite to the area with the test pads. Do not open the tube containing the test strips unless you are ready to conduct the test. Remove only as many test strips as required and reseal the tube immediately with the cap. The cap contains a desiccant. Perform testing in ambient temperature (20±5°C).

Collection of Urine Sample:

- Each urine sample must be collected in a clean container made of polypropylene. Only freshly collected urine must be used. If urine samples are collected by urinary catheter, the maximum time between collection in the catheter and testing is two hours.
- Shake the container lightly before dipping the test strip.
- Use urine within two hours after collection.

Assay Procedure - Use of Test Strips:

- Immerse the test strip for 2 to 3 seconds into the urine sample so that all test pads are completely covered by the urine.
- After removing the test strip from the urine, some liquid may be attached at the borders or edges of the test strips. Wipe off the excess urine on a tissue to absorb runoff. The test pads should not be touched.
- Place the test strip on a flat surface so that you can see the test pads, and wait for 10 minutes (incubation time of tests). Use a timer to check the time.
- After 10 minutes incubation time immediately compare the test pads by naked eye to the corresponding colour scales on the label of the tube container. Refer to the next section regarding the visual determination of the colours.

Visual Determination and Interpretation of Colours of DOAC Dipstick

Colours must be compared to colours of the colour scale printed on the label of the container by naked eye.

Test Pad 1 (Creatinine):

- The colour of pad 1 corresponds to colours "norm." on the tube label → creatinine in urine is normal. Pad 3 and pad 4 can be evaluated.
- The colour of pad 1 is darker than the colour "norm." on the tube label → high creatinine does not affect DOAC excretion into urine. Pad 3 and pad 4 can be evaluated.
- The colour of pad 1 is "low" or lighter than the respective colour of the tube label → creatinine in urine is low, indicating renal insufficiency. **Colours of pad 3 and pad 4 may be false negative.**

Test Pad 2 (Urine colour):

- The colour of the pad is as white as the respective colour marked "norm." on the tube label → the results of pad 1, 3 and 4 are valid.
- The colour of the pad is darker than the colour printed on the tube label → colours of pad 1, pad 3 and pad 4 may be distorted. The test is invalid.

Test Pad 3 (Medication Apixaban, Edoxaban, Rivaroxaban):

- The colour of pad 3 clearly is yellow as the respective colour marked "neg." on the tube label → direct oral Factor Xa inhibitor is absent in the urine sample.
- The colour of pad 3 is less yellow than the colour marked "neg." on the tube label, thus the result is "pos." → direct oral Factor Xa inhibitor is present in urine.
- The colour of pad 3 is white as the respective colour marked "pos." on the tube label → direct oral Factor Xa inhibitor is present in urine.

Test Pad 4 (Medication Dabigatran):

- The colour of pad 4 is ochre as the respective colour marked "neg." on the tube label → direct oral Thrombin inhibitor is absent in the urine sample.
- The colour of pad 4 is between the ochre colour marked "neg." and the rose colour marked "pos." on the tube label → direct oral Thrombin inhibitor is present in urine.
- The colour of pad 4 is rose as the respective colour marked "pos." on the tube label → direct oral Thrombin inhibitor is present in urine.

If pad 3 and pad 4 are both "pos.", the test is invalid because it is unlikely that a person is treated with both types of DOACs.

Limitations

The test results are qualitative. No quantitative interpretation should be made based on the test results. The results always must be interpreted and evaluated in connection with other clinical information by a physician before diagnosing. No treatment decisions should be made solely based on the outcome of a DOAC Dipstick analysis. Additional laboratory analysis (e.g. determination of blood coagulation parameters) may be required. A recently voided urinary bladder may impact the detectable DOACs in the urine. Persons with colour vision deficiency or colour blindness must not perform the DOAC Dipstick test.

Expected Values and Reference Ranges

Creatinine – Reference Range: 0.25 – 3.0 g/l, (2.2 – 26.5 mmol/l) (Ref.: Needleman).

DOACs – Normal values are below 5 ng/ml (LC-MS/MS-method). Patients under DOAC treatment typically display values above 200 ng/ml (Ref.: Schreiner).

Performance Characteristics

The visual analysis of the colours of the DOAC Dipstick corresponds to results "negative" and "positive". The cut-off value of Apixaban, Edoxaban and Rivaroxaban for "negative" is < 100 ng/ml, and the cut-off value for "positive" is > 200 ng/ml. The cut-off value of Dabigatran for "negative" is < 50 ng/ml and for "positive" >125 ng/ml. In the ranges between the cut-off values the colours of the results for the DOACs may be identified as either "negative" or "positive". Data were obtained in artificial urine and normal human urine spiked with concentrations of DOACs ranging from 0 to 1500 ng/ml. Using urine samples of patients treated with Apixaban, Edoxaban, Rivaroxaban, and Dabigatran, sensitivities and specificities were >95% (Ref.: Harenberg 2020). Concentrations of DOACs in urine are higher due to the lower volume of urine compared to the volume of distribution of DOACs in blood, and were below 5 ng/ml (LC-MS/MS-method) in patients not treated with DOACs and typically above 200 ng/ml for all DOACs during treatment and 12 or 24 hours after intake of medication (mean values 5.600 ng/ml Dabigatran, 2.700 ng/ml Rivaroxaban, 1.800 ng/ml Apixaban, n=29 each) (Ref.: Schreiner).

Interferences

No information is reported in literature on drug-drug and drug-other compounds interactions, except the coloured compounds in urine as described above. The number of interactions tends to be low to very low due to the high specificity of the enzymes with the respective substrates. No interaction occurs between the components of pad 3 with pad 4 and vice versa. Heparins do not react on pads 3 and 4 (Ref.: Harenberg 2017). Coloured compounds in urine such as bilirubin, urobilinogen and blood (macrohaematuria) may modify the results of pad 1, pad 3 and pad 4. The impact of the colour of the urine sample can be assessed by pad 2 (see above for interpretation). Creatinine: In urines with high buffering capacity false negative reading may be obtained. With high concentrations of ketone bodies (>50 mmol/l) false positive reading may occur. Blood >2000 Ery/µl may cause false positive results.

Storage Conditions

The recommended storage temperature of the closed test-strip tube container is between +2 °C and +30 °C. DOAC Dipsticks may be used until the expiry date printed on the label. After opening, the tube must be closed tightly immediately after removing the strips used for evaluation, and the closed tube has to be stored away from direct sunlight at temperatures not exceeding +30 °C. Storage above +30 °C will adversely affect the stability and test performance of the product. After first opening the tube, the DOAC Dipsticks remaining in the tube can be used for a maximum of 3 month. Do not use any DOAC Dipsticks remaining in the tube after 3 months after first opening.

References

- CLSI. Urinalysis; Approved Guideline - Third Edition. CLSI document GP16-A3; 2009.
Schreiner R et al. Res Pract Thromb Haemost 2017; 1(Suppl.1): PB 491.
Harenberg J et al. Thromb Haemost 2020, 120:132-40.
Harenberg J et al. Res Pract Thromb Haemost 2017; 1(Suppl.1): PB 454.
Harenberg J et al. Clin Chem Lab Med 2016; 54: 275-83.
Du S et al. Clin Chem Lab Med 2015; 53: 1237-47.
Harenberg J et al. Semin Thromb Hemost 2015; 41: 228-36.
Favaloro EJ et al. Semin Thromb Hemost 2015; 41: 208-27.
Harenberg J et al. Thromb J 2013 Aug 1; 11(1): 15.
Needleman SB et al. J Forensic Sci 1992; 37: 1125-33.

W17.5-3-AU-Rev01 – Date of revision: 2020-05-25

SYMBOLS USED

- CE Mark - Device complies with the Directive 98/79/EC
- LOT Lot number
- REF Catalogue number
- IVD In vitro diagnostics
- Manufacturer
- Storage temperature
- Expiry date
- Content is sufficient for <n> tests
- Follow carefully instructions for use

DOASENSE™ GmbH, Waldhofer Str. 102, 69123 Heidelberg, Germany
Email: info@doasense.de, www.doasense.de