

New programs for 2022 catalogue

## DOAC Urine Dipstick Program

Program suitable for detection of Direct Oral Anticoagulants in urine samples using (DOASENSE) urine dipstick.

**Catalogue number:** 11540102

**Frequency (No of surveys/year):** 2

**Cases/samples per survey:** 2

**Schedule:** May, November

**Sample type:**

Lyophilised urine. Reconstitution consumables provided

**Measurands:**

Creatinine, Thrombin Inhibitor, Xa Inhibitor

**Assessment:**

Qualitative interpretation of measurands

**Storage and handling:**

Refrigerate samples (2-8°C) upon receipt

**Transport:**

Exempt human, shipped and stable at ambient temperatures

**Shipment frequency:**

Packaged kit box containing all samples and consumables for full schedule of testing. Dispatched once a year, before first survey open.

### DOAC Urine Dipstick

- Rapid point-of-care qualitative strip test for detecting Direct Oral Anticoagulants (DOAC Dipstick) in urine.
- Excludes presence of clinically relevant DOAC (both factor Xa inhibitors and thrombin inhibitors) concentrations.
- Useful when medication history is unavailable and where laboratory testing for DOACs is not available.
- The DOAC Dipstick test can provide rapid decision support for patients with severe bleeding or patients with indication for an intervention, e.g., in trauma surgery or neurology.

